
**BW Series Compact Bench Scale**

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**Balance de table compacte, série BW**

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**Bilancia da banco Compact serie BW**

**Ohaus Corporation, 19A Chapin Road, PO Box 2033, Pine Brook, NJ 07058 USA**

**Declaration of Conformity** We, Ohaus Corporation, declare under our sole responsibility, that the scale models listed below marked with "CE" are in conformity with the directives and standards mentioned.

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Ted Xia  
 President  
 Ohaus Corporation  
 Pine Brook, NJ USA

Urs Müller  
 General Manager  
 Ohaus Europe  
 Grenensee, Switzerland

## Important notice for verified weighing instruments

Weighing Instruments verified at the place of manufacture bear one of the preceding mark on the packing label and the green 'M' (metrology) sticker on the descriptive plate. They may be put into service immediately.

Weighing Instruments to be verified in two stages have no green 'M' (metrology) on the descriptive plate and bear one of the preceding identification mark on the packing label. The second stage of the initial verification must be carried out by the approved service organization of the authorized representative within the EC or by the national weight & measures (W+M) authorities.

The first stage of the initial verification has been carried out at the manufacturers work. It comprises all tests according to the adopted European standard EN 45501:1992, paragraph 8.2.2.

If national regulations limit the validity period of the verification, the user of the weighing instrument must strictly observe the re-verification period and inform the respective W+M authorities.

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In conformance with the European Directive 2002/96/ EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.

If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environmental protection.

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### **ISO 9001 Registration**

In 1994, Ohaus Corporation, USA, was awarded a certificate of registration to ISO 9001 by Bureau Veritus Quality International (BVQI), confirming that the Ohaus quality management system is compliant with the ISO 9001 standard's requirements. On May 15, 2003, Ohaus Corporation, USA, was re-registered to the ISO 9001:2000 standard.

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# 1. INTRODUCTION

This manual contains installation, operation and maintenance instructions for the BW Series Compact Bench Scale. Please read the manual completely before installation and operation.

## 1.1 Product Description

The BW Series Compact Bench Scale is an economical portioning, packing and general weighing scale designed to provide fast and accurate weighing performance in washdown environments. With an internal sealed rechargeable lead-acid battery, IP65 housing and high-speed display update, the BW Series Compact Bench Scale offers portability, durability and productivity in one practical design.

## 1.2 General Features

- Capacity models: 1.5kg, 3kg, 6kg, 15kg
- Maximum displayed resolutions from 1:6,000-7,500
- Type approved resolution 1:3,000
- High-profile ABS housing, stainless steel pan (pan is dishwasher-safe)
- IP65 dust and water jet protection
- 4-adjustable leveling feet with integrated front-panel leveling indicator
- Bright red LED weight display: 6-digit, 7-segment, 14 mm high characters
- 4-button control panel: On/Off, Check, Tare, Zero
- 2-weighing units: kg, g
- Battery charge indication
- AC power or 15-hour rechargeable lead-acid battery
- Over-Under checkweighing function with audible alarm

## 1.3 Safety Precautions

For safe and dependable operation of this scale, please comply with the following safety precautions:

- Verify that the input voltage printed on the AC Adapter and the plug type matches the local AC power supply.
- Make sure that the power cord does not pose an obstacle or tripping hazard.
- Disconnect the scale from the power supply when cleaning the scale.
- Do not operate the scale in hazardous or unstable environments.
- Do not drop loads on the platform.
- Use only approved accessories and peripherals, as available.
- Operate the scale only under ambient conditions specified in these instructions.
- Service should be performed by authorized personnel only.

## 2. INSTALLATION

### 2.1 Unpacking

Unpack and verify that the following components have been included:

- BW Scale Unit
- Weighing Pan
- AC-DC Adapter
- Instruction Manual


Save the packaging material. This packaging ensures the ideal protection for the storage or transport of the product.

### 2.2 Transport Protection

#### BEFORE USING THE SCALE:

Remove the 3 shipping protection screws on the spider sub-platform, then remove the 3 plastic column inserts under the spider. Re-install the protective components if shipping the scale to ensure proper transport protection.


### 2.3 Selecting the Location

Operate the scale on a firm, level surface. Avoid locations with rapid temperature changes, excessive air currents, vibrations, electromagnetic fields, heat or direct sunlight.

### 2.4 Leveling the Scale

Adjust the leveling feet\* until the bubble is centered in the circle of the level indicator (located on the front panel). Ensure that the scale is level each time its location is changed.

\* **IMPORTANT:** To allow proper air flow and minimize condensation inside the scale, unscrew the 2 rear leveling feet slightly, then level the scale accordingly for normal operation. However, when cleaning the scale, all 4 leveling feet must be screwed in tightly to avoid entry of water inside the scale.

## 2.5 Connecting Power

### 2.5.1 AC Power

Verify that the intended AC power source matches the AC adapter rating. Unplug\* the rubber access cap underneath the scale, connect the adapter to the power receptacle through the cut-out, then plug the AC adapter into a power outlet.

\* **NOTE:** When cleaning the scale, turn off the unit, remove the AC adapter and plug the rubber access cap securely to avoid entry of water inside the scale.

### 2.5.2 Battery Power

The scale will begin charging with the AC adapter connected accordingly, and automatically switches to battery operation if there is a loss of power. Normal operating time on a full charge is 50 hours on a 15-hour charge cycle. DO NOT continuously charge the scale longer than 24 hours.

Low battery charge is indicated by the flashing LED over the  indicator. AC power or a charged battery is indicated by the illuminated LED over the  indicator.

#### NOTES:

- To maximize battery power, the scale can be set up to enter Sleep mode (only decimal points are displayed) after 2 minutes of inactivity, or automatically turn off after 15 minutes of inactivity (Section 5.2).
- Charging the scale must be performed in a dry environment.
- Plug the rubber access cap to the adapter receptacle to ensure the scale is protected when operated in a wet environment.
- If the scale is not used within 3 months, recharge and discharge the battery fully 3 times to maximize the battery condition.
- Dispose of the lead-acid battery according to local laws and regulations.

#### CAUTION



**BATTERY IS TO BE REPLACED ONLY BY AN AUTHORIZED OHAUS SERVICE DEALER. RISK OF EXPLOSION CAN OCCUR IF REPLACED WITH THE WRONG TYPE OR CONNECTED IMPROPERLY.**

## 2.6 Initial Calibration

When the scale is operated for the first time, a Span Calibration is recommended to ensure accurate weighing results. Before performing the calibration, be sure to have the appropriate calibration weights. Refer to Section 6 for Span Calibration procedures.

### 3. OVERVIEW OF CONTROLS AND FUNCTIONS

- 
  - Turns the scale on
  - Long-press turns the scale off
  - In SETUP Mode, “Yes” accepts/confirms the displayed setting
  - In CHECK Mode, accepts/confirms the displayed value
  
- 
  - Inputs the weight of the object on the weighing pan as a Tare value
  - In SETUP Mode, “No” toggles through the available settings
  - In CHECK Mode, increments the value of the selected digit
  
- 
  - Zeros the display
  - With the scale off, long-press together with **ON/OFF** enters User Setup Mode
  - In SETUP Mode, “Exit” ends setup and exits to “Save”
  - In CHECK Mode, ends Over/Under setup and exits to weighing mode
  
- 
  - Initiates CHECK Over/Under setup
  - In SETUP Mode, “Back” reverts to the previous parameter
  - In CHECK Mode, shifts to the next digit
  
- 
  - LED annunciator indicates Center of Zero
- 
  - LED annunciator indicates power is AC or battery is fully charged
- 
  - LED annunciator indicates displayed weight is less than the preset Under value
- 
  - LED annunciator indicates displayed weight is more than the preset Over value
- 
  - LED annunciator blinks to indicate the battery voltage is low
- 
  - LED annunciator indicates Net (Gross minus Tare) value
- 
  - LED annunciator indicates the active weighing unit

## **4. OPERATIONS**

### **4.1 Turning the Unit On and Off**

With the scale off, press **ON/OFF** to turn on. Make sure the weighing pan is empty before turning the scale on. Allow 15-30 minutes for the scale to warm up before use. With the scale on, long-press **ON/OFF** to turn off.

### **4.2 Manual Tare**

Place a container on the weighing pan, then press **TARE**. The display will show a net value of "0" as the Tare weight is registered. Remove the container, and the resulting negative value displayed will be equivalent to the Tare weight. To clear the Tare value, press **TARE** with the pan empty.

### **4.3 Zero Operation**

Press **ZERO** to zero the stable weight display (the weight must be within  $\pm 2\%$  of scale capacity). Center of Zero is indicated by the LED annunciator.

### **4.4 Checkweigh Target Setting**

Compares the weight of an item against a pre-set Over-Under range.

To activate Checkweighing Mode, refer to user setup (Section 5.2).

- Press **CHECK** to initiate entering Over-Under parameters.
- The value of the "OVER" target limit will be displayed first, with the active digit flashing (the "OVER" LED annunciator is lit).
- Enter the desired "OVER" limit value by pressing **TARE 0 9** to increment the active digit or **CHECK** to shift to the next digit.
- Press **ON/OFF** to accept the displayed value and proceed to the "UNDER" target limit (the "UNDER" LED annunciator is lit).
- Enter the desired "UNDER" limit value by pressing **TARE 0 9** to increment the active digit or **CHECK** to shift to the next digit.
- Press **ZERO** to accept the displayed value and exit into checkweighing mode.

#### **NOTES:**

- "OVER" or "UNDER" LED annunciators will light if the weight of an item placed on the pan is over or under the pre-set limits, respectively. "OVER" and "UNDER" LED annunciators will light simultaneously if the weight is within the pre-set limits.

- The LED indication is complemented by an audible alert feature (Section 5.2).
- Over-Under values are retained when the unit is turned off.
- Over-Under values are not retained when the weighing unit or resolution setup are changed.

## 5. SCALE SETTINGS

### 5.1 Accessing the User Setup Mode

With the scale off, press **ON/OFF** and **ZERO** simultaneously. Once the scale turns on, release **ON/OFF** while continuously pressing **ZERO** until “**SETUP**” is displayed. To navigate in the User Setup Mode:

- Press **ON/OFF-Yes** to accept the displayed setting
- Press **TARE-N'o** to toggle through available settings
- Press **ZERO-Exit** to end setup and proceed to “Save”
- Press **CHECK-Back** to revert to the previous parameter

### 5.2 User Setup Parameters

Parameters	Description	Settings
S1	Scale automatically turns off after 15 minutes of inactivity.	<b>ON</b> - enables this function. <b>OFF</b> -disables this function.
S2	Scale enters sleep mode after 2 minutes of inactivity.	<b>ON</b> - enables this function. <b>OFF</b> -disables this function.
S3	CHECK activation and Alert Range	<b>0</b> - Check mode is off.  <i>Scale beeps when displayed weight is:</i> 1-...within the preset Over and Under values 2-...outside the preset Over and Under range 3-...within the preset Over and Under values (Over-Under LED's do not light within preset range) 4-...outside the preset Over and Under range (Over-Under LED's do not light within preset range) 5- Alarm is off (Over-Under LED's do not light within preset range) 6- Alarm is off
S4	Averaging Level	0, 1, <b>2</b> , 3 (most filtering, longer stabilization time) 4, 5, 6, 7 (corresponds to 0, 1, 2, 3 but with slower display update)
S5	Weighing Unit	<b>kg</b> , g
S6	LED Brightness Level	<b>0</b> , 1, 2, 3 (dimmmst).
S7	Display Type	<b>OFF</b> - displays weight continuously. <b>ON</b> - displays first stable weight.
S8	Resolution	<b>OFF</b> - 3000 scale divisions <b>ON</b> - 6000/7500 scale divisions
Save	Save or Abort	<b>Save</b> - saves changes and reverts to weighing mode. Abort - reverts to weighing mode without saving changes.

**Note:** Bold settings are defaults.


## 6. CALIBRATION AND SEALING

### 6.1 Calibration

For best results, calibrate the scale at regular intervals. Temperature changes, geographic gravity variations, altitude changes and abuse are few reasons why a scale may need recalibration.

**Note:** The following procedure requires the use of hand tools and calibration weights. If unable to perform the calibration, contact an authorized Ohaus service organization.

1. Remove the battery cover to access the two calibration pins located in a slot at the bottom of the scale, then short the calibration pins (a small slotted screwdriver may be used).

**Note:** When accessing the bottom of the scale, avoid placing the scale upside down on the pan or sub-platform. Place the scale on its side.

2. **SETUP** then the first parameter **CIY** is displayed.
3. Press **ON/OFF-Yes** successively to move through the succeeding parameters: "dEF NO", "S1", "S2", "S3", "S4", "S5", "S6", "dSP", "S8", "S9", "S10", "S11", "S12", "bAt", "GEO".

**Note:** Do not change the settings of these parameters. These are service setup parameters intended for authorized servicing use only.

4. After "GEO", **CAL NO** is displayed. Press **TARE-No** to switch the setting to **CAL YES**.
5. **CAL YES** is displayed. Press **ON/OFF-Yes** to enter calibration sequence.

6. "kg" is displayed. This is the Calibration Weight Unit. Press **ON/OFF-Yes** to accept and move on to the next parameter.
7. "xxxx" is displayed. This is the Full Scale Capacity of the scale, where: "1500" = 1.5kg, "3000" = 3kg, "6000" = 6kg, "15000" = 15kg. Make sure the displayed value matches the scale model capacity then press **ON/OFF-Yes** to accept and move on to the next parameter.
8. "- - - - -" is displayed. Make sure the pan is empty and stable at this time then press **ON/OFF-Yes**. The display will count down from 5 to 1 as the scale zero is registered.
9. "xxxx" is displayed. This is the Calibration Point corresponding to 2/3 of the Full Scale Capacity value. Place the indicated calibration weight on the center of the pan. With the scale stable, press **ON/OFF-Yes**. The display will count down from 5 to 1 as the scale registers the calibration weight.
10. "xxxx" is displayed. This is the Calibration Point corresponding to the Full Scale Capacity value. Place the indicated calibration weight on the center of the pan. With the scale stable, press **ON/OFF-Yes**. The display will count down from 5 to 1 as the scale registers the calibration weight.
11. "SAVE" is displayed when calibration was registered. Remove the calibration weight, then press **ON/OFF-Yes**. The scale will display a self-test then return to normal weighing mode.


**Note:** At any time during the calibration process, press **CHECK-Back** to revert to the previous parameter, or long-press **ON/OFF** to turn off the scale and abort the process without saving any changes.

## **6.2 Type Approval Sealing**

After the scale has been found to comply with local regulations by relevant approving authorities, the scale may be sealed in the following manner:


### **Lead and Wire Seal.**

Secure sealing plate with through-hole screws. Run the wire seal into the through-holes then press the lead seal in place.

### **Sticker Seal.**

Secure sealing plate with regular pan-head screws then apply the sticker over one of the screws.

## **7. CLEANING AND MAINTENANCE**

Comply with safety precautions when cleaning and maintaining the scale

- Turn off and disconnect the scale from the power supply during cleaning.
- **DO NOT** open the scale; no user repairable components inside. Service should be performed by authorized personnel only.
- **DO NOT** place the scale upside down. If accessing the bottom, place the scale on its side.
- **DO NOT** carry the scale by the pan or sub-platform. Use the handholds on the sides of the scale housing.
- **DO NOT** allow liquids to enter the scale.
- **DO NOT** use high pressure when cleaning the scale.
- Inspect the scale regularly. Ensure that the housing, pan and sub-platform are clean and free from foreign material.
- Wipe the exterior components of the scale using a soft cloth dampened with water and a mild detergent (**DO NOT** use acids, alkalis, strong solvents or abrasive cleaning agents).
- Screw in all 4 leveling feet tightly to avoid entry of water inside the scale during cleaning. However, when the scale is put back into operation, remember to loosen the two rear feet (and adjust leveling accordingly) to allow proper air flow and minimize condensation inside the scale.
- Plug the rubber cap to the adapter receptacle underneath the scale to avoid entry of water inside the scale.
- The pan is dishwasher safe and may be placed with normal loads. All other parts must not be placed in a dishwasher.


## 8. TROUBLESHOOTING

The following table lists common problems, possible causes and remedies. If the problem persists, contact Ohaus or your authorized Ohaus dealer.

<b>Symptom</b>	<b>Possible Causes</b>	<b>Remedy</b>
Scale will not turn on	AC power not connected to scale or to outlet. No power in outlet. Battery discharged.	Connect scale to power. Connect scale to power and charge battery. Connect to different power outlet. Replace battery.
Battery fails to charge fully.	Battery is defective or past its useful life.	Have the battery replaced by an authorized Ohaus service dealer.
Weight reading does not stabilize.	Unstable environment. Interference under pan or sub-platform.	Ensure a stable environment. Ensure that the weighing pan and sub-platform are unobstructed.
Scale does not display accurately.	Improper calibration	Calibrate the scale using proper calibration weights.
“nnnnn” displayed	Overload condition: Load on the pan exceeds scale capacity +9d	Reduce the load on the pan.
“uuuuu” displayed	Underload condition: No load on the pan, exceeds scale zero -9d	Install the pan on the sub-platform. Re-zero the scale.
“- - - - -” displayed	Unstable condition. Cannot establish zero.	Situate in a stable environment. Relieve internal pressure from the scale by loosening the leveling feet. Remove the weight on the pan that may exceed the power-on zero range.
“E 11” displayed	RAM error	Turn the scale off and on again.
“E 16” displayed	ROM error	Recalibrate the scale.
“E 18” displayed	EEPROM error	Return scale for service
“E 48” displayed	Checkweigh target setting error	Review setup and re-enter correct Over-Under values
“Error” displayed	Software running error	Turn the scale off and on again. Replace battery.

## **9. TECHNICAL DATA**

### **9.1 Drawings**

## 9.2 Specifications

Model	BW1.5	BW3	BW6	BW15
Capacity x Readability	1.5 kg x 0.2 g	3 kg x 0.5 g	6 kg x 1 g	15 kg x 2 g
Maximum Displayed Resolution	1:7,500	1:6,000	1:6,000	1:7,500
Type Approval Resolution *	1:3,000			
Load Cell Capacity	3.5kg	6.5kg	15kg	30kg
Construction	ABS plastic housing with stainless steel pan			
Protection	IP65			
Pan Dimensions	230W x 180D mm			
Leveling Components	4 x adjustable feet with integral panel mounted level bubble indicator			
Display	1 x red LED; 6-digit, 7-segment, 14H mm characters			
Weighing Units	kg, g			
Keyboard	4 x tactile mechanical keys: ON/OFF, CHECK, TARE, ZERO			
Calibration Weight	2/3 and Full Scale Capacity			
Over Range Capacity	Full Scale Capacity +9d			
Tare Range	Full Capacity by Subtraction			
Zero Range	Power-on: $\pm 10\%$ / Keyboard: $\pm 2\%$ Full Scale Capacity			
Stabilization Time	within 2 seconds			
Power	9VDC 500mA AC Adapter / 6V 5Ah rechargeable lead acid battery			
Battery Life	50 hours continuous use with 15 hour recharge time			
Safe Overload	150% Full Scale Capacity			
Operating Temperature Range	-10 ° C to +40 ° C, 10% to 85% humidity, non-condensing			
Storage Temperature Range	-25 ° C to +50 ° C, 10% to 85% humidity, non-condensing			
Scale Dimensions	230W x 310D x 127H mm			
Shipping Dimensions	295W x 381D x 190H mm			
Scale Unit Weight	4.4 kg			
Shipping Weight	5.7 kg			

\* OIML 3000e: R76/1992-NL 1-03.49

## **LIMITED WARRANTY**

Ohaus products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period, Ohaus will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to Ohaus.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than Ohaus. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by Ohaus Corporation. Ohaus Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact Ohaus or your local Ohaus dealer for further details.

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P/N 80251077 A © Ohaus Corporation 2006, all rights reserved / todos los derechos reservados / tous droits réservés / Alle Rechte vorbehalten / tutti i diritti riservati.